






# ORAL FLUID INITIAL AND CONFIRMATORY TEST REQUIREMENTS

With the addition of lab-based oral fluid to federal workplace drug testing comes a new era for the industry. Although SAMHSA has aligned as much of the oral fluid guidelines as possible with the pre-existing urine guidelines, some things must change, such as initial and confirmatory test requirements.

*The below information details initial and confirmatory test requirements for lab-based oral fluid only.*

	Initial Test Analyte	Initial Test Cutoff	Confirmatory Test Analyte	Confirmatory Test Cutoff Concentration
	Marijuana (THC)	4 ng/mL	THC	2 ng/mL
	Cocaine/ Benzoyllecgonine	15 ng/mL	Cocaine Benzoyllecgonine	8 ng/mL
	Codeine/ Morphine	30 ng/mL	Codeine Morphine	15 ng/mL
	Hydrocodone/ Hydromorphone	30 ng/mL	Hydrocodone Hydromorphone	15 ng/mL
	Oxycodone/ Oxymorphone	30 ng/mL	Oxycodone Oxymorphone	15 ng/mL
	6-Acetylmorphine	4 ng/mL	6-Acetylmorphine	2 ng/mL
	Phencyclidine	10 ng/mL	Phencyclidine	10 ng/mL
	Amphetamine/ Methamphetamine	50 ng/mL	Amphetamine Methamphetamine	25 ng/mL
	MDMA/MDA	50 ng/mL	MDMA MDA	25 ng/mL

*This information is provided for educational purposes only. Reader retains full responsibility for the use of the information contained herein.*

©2010-2019 The Current Consulting Group, LLC – No portion of this article may be reproduced, retransmitted, posted on a website, or used in any manner without the written consent of the Current Consulting Group, LLC. When permission is granted to reproduce this article in any way, full attribution to the author and copyright holder are required.